



¿Qué es el material genético,
cómo obtenerlo y para qué?

Nataly Ruiz Quiñones

nruizq@gmail.com



Breve historia de la molécula de DNA

1869- Friedrich Miescher.

- Primera extracción de DNA.
- Nucleina (Miescher et al, 1971)



<https://www.advancedsciencenews.com/wp-content/uploads/2019/04/Capture-2-233x300.png>

1879, 1890- Albrecht Kossel
(premio Nobel 1910)

- Identificó las bases nitrogenadas
- Presencia de azúcar.




<https://iiif.wellcomecollection.org/image/V0026654.jpg/full/800%2C/0/default.jpg>

1889- Richard Attmann

- Ácido nucleico




<https://s3.amazonaws.com/s3.amazonaws.com/public/uploads/photos/5356276/richard-altmann-1852-1900.jpg>



Breve historia de la molécula de DNA

1929- Phoebus Levene.

- Ribosa, desoxirribosa y fosfato.
- Nucleotidos



<https://images.fineartamerica.com/image-s-medium-large-5/phoebus-levene-national-library-of-medicine.jpg>

1950- Erwin Chargaff

- Nucleotidos inter e intraespecífico.
- $A=T$, $G=C$




<https://media.sciencephoto.com/image/c0373020/800wm>

1952- Alfred Hershey y Martha Chase

- Carácter hereditario del DNA




https://content.dnalc.org/content/c16/16406/16406_hershey_chase.jpg



Breve historia de la molécula de DNA

1953- Rosalind Franklin y Maurice Wilkins


- Identifican la doble hélice del DNA




<https://www.kcl.ac.uk/ImportedImages/0Prospectus/kings-in-time/dna-pic-2.jpg>

1953- James Watson y Francis Crick




- Dilucidan la estructura del DNA




https://www.sciencehistory.org/sites/default/files/styles/rte_full_width/public/watson-crick-dna-model.jpg?itok=Qa7645Jc



Conformación del DNA




https://pm1.narvii.com/6594/da2cf51b0658c0e5cbe2d78877dca3853dfe54f_hq.jpg



(a) Molécula de DNA lineal de doble cadena

Molécula de DNA circular

(a)




(b)


Fuente: Trudy McKee, James R. McKee: *Bioquímica. Las bases moleculares de la vida*, 5e: www.accessmedicina.com
Derechos © McGraw-Hill Education. Derechos Reservados.

<https://www.google.com/url?sa=i&url=https%3A%2F%2Faccessmedicina.mhmedical.com%2FContent.aspx%3Fbookid%3D1960%26sectionid%3D148097228&psig=AOvVaw19inbtzRxRxbF12jJlezZt&ust=1576204616824000&source=images&cd=vfe&ved=0CAIQjRxqFwoTCIDj9piPr-YCFQAAAAAdAAAAABAD>



Conformação del RNA



<https://player.slideplayer.com.br/26/8862880/data/images/img11.jpg>





Conformación del DNA y el RNA






Dogma central de la biología molecular



1957: Francis Crick



Procesos de obtención del material genético



¿Cómo seleccionar el mejor protocolo de extracción?



¿Cómo seleccionar el mejor protocolo de extracción?

4. Finalidad del material

Pureza, calidad y cantidad

Integridad (Téc. molecular)



5. Método a utilizar

Disponibilidad económica


Capacitación del personal





6. Preservación del material genético


Uso inmediato

A largo plazo



Procesos de obtención del material genético





	Origen de muestra								
	CC	TM	Sa	Ba	Fu	Al	Le	Pl	In
microbios									
Bacterial Genomic DNA Mini-prep Kit (BayGene)				*					
QIAprep Spin Miniprep Kit (Qiagen)				*					
QIASymphony Virus/Bacteria Kits (Qiagen)				*					
ZR Fecal DNA Mini Prep (Zymo Research)			*	*					
Plasmid Maxi Kit (Qiagen)				*					
BACMAX DNA purification kit (Epicentre Biotechnologies)				*					
PowerMax Soil DNA Isolation Kit (MO BIO Laboratories)				*	*	*			
PowerSoil DNA isolation kit (MO BIO Laboratories)				*	*	*			
Células y tejidos de mamíferos									
AccuPrep Genomic DNA Extraction Kit (Bioneer)	*	*	*						
Arcturus DNA Extraction Kit (Arcturus)	*	*							
GFX Genomic Blood DNA Purification Kit (GE Healthcare)	*		*						
DNA Isolation Kit for mammalian blood (Roche)			*						
InnuPrep DNA minikit (AJ Innuscreen)	*	*							
QIAamp DNA mini kit (Qiagen)	*	*	*						
AllPrep DNA/RNA Mini Kit (Qiagen)	*	*							
Agencourt DNAdvance Kit (Beckman Coulter)	*	*							

	Origen de muestra								
	CC	TM	Sa	Ba	Fu	Al	Le	Pl	In
plantas									
NucleoSpin 8 Plant and NucleoSpin 96 Plant II, Clontech								*	
DNeasy 96 Plant Kit (Qiagen)					*			*	
Nucleon PhytoPure Genomic DNA Extraction Kits (GE Healthcare)					*			*	
Células de mamíferos y microbios									
DNA Isolation Kit for cells and tissues (Roche)	*	*		*			*		
Purelink Genomic DNA extraction kit (Invitrogen)	*	*	*	*					
DNeasy Blood and Tissue Kit (Qiagen)	*	*	*	*			*		*
Genomic DNA from Tissue kit (Macherey Nagel)	*	*	*	*			*		
GeneJET Genomic DNA Purification Kit	*	*	*	*			*		
FastDNA SPIN Kit (MP Biomedicals)	*	*		*	*	*	*	*	*
Células de mamíferos, microbios y plantas									
ArchivePure DNA purification kit (5Prime)	*	*	*	*			*	*	
DNA Isolation Kits (BioBasic)	*	*	*	*	*		*	*	
FastDNA Kit (MP Biomedicals LLC)	*	*	*	*	*	*	*	*	*
DNAzol® Reagent (Invitrogen)	*	*	*	*	*	*	*	*	*
Easy-DNA® Kit (Invitrogen)	*	*	*	*			*	*	
Wizard® Genomic DNA Purification Kit (Promega)	*	*	*	*			*	*	*
sangre									
DNA Isolation Kit for mammalian blood (Roche)			*						
InnuPREP Blood DNA Mini Kit (AJ Innuscreen)			*						

Tabla 1. Diferentes kits de extracción de DNA y RNA. Tomado de: <http://www.labome.es/method/DNA-Extraction-and-Purification.html>

CC: células en cultivo; TM: tejido de mamíferos; Sa: sangre; Ba: bacterias; Fu: hongos (fungi); Al: algas; Le: levaduras; Pl: plantas; In: insectos

<i>Productos de PCR</i>	<i>Extracción de gel</i>
Wizard® PCR Preps DNA Purification System (Promega)	MinElute Gel Extraction Kit (Qiagen)
QIAquick PCR Purification Kit (Qiagen)	Zymoclean™ Gel DNA Recovery Kit (Zymo Research)
MinElute PCR Purification Kit (Qiagen)	<i>Otros Kits</i>
GenElute™ PCR Clean-Up (Sigma-Aldrich)	NEBNext DNA Sample Prep Reagent Set 1 (New England Biolabs)
PureLink® PCR Purification Kit (Life Technologies)	GS FLX Titanium Rapid Library Preparation Kit (Roche)
GeneJet PCR Purification Kit (Thermo Scientific)	

Tabla 2. Otros kits de purificación de DNA . Tomado de: <http://www.labome.es/method/DNA-Extraction-and-Purification.html>

How to isolate total RNA using magnetic bead-based technology

Applied Biosystems™ MagMAX™ mirVana™
Total RNA Isolation Kit
High-throughput, pure, concentrated RNA

https://www.youtube.com/watch?v=h6NvCd94eng&list=PLGlVFEwL2wDFs3L_mbHpyfNbRrj5W3A0Q&index=4

How to isolate RNA from tissue or cells

Invitrogen™ Ambion™ TRIzol™ Plus RNA Purification Kit

Ultrapure total RNA within an hour,
even from difficult samples. Ideal for
RNA sequencing, RT-PCR, and microarrays.

https://www.youtube.com/watch?v=L5XIB8KkKt4&list=PLGlVFEwL2wDFs3L_mbHpyfNbRrj5W3A0Q&index=8

How to purify microbial and host DNA from stool samples

Invitrogen™ PureLink™ Microbiome DNA Purification Kit
Fast, high-quality DNA isolation

https://www.youtube.com/watch?v=7muQSoATQV8&list=PLGlVFEwL2wDFs3L_mbHpyfNbRrj5W3A0Q&index=6

